



Disposable Surgical Gown SMMS, Non-Sterile

  **MADE
IN TURKEY**

Product Details

Product Description : SMMS Surgical Gown

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSMG-SMMS08

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21050501

Product Name: Disposable non-sterile surgical gown

Type : SMMS Surgical Gown

Usage: Standard surgical gown

Sterilization: Non-sterile

Weight: 35-45 g / m²

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SMMS Surgical Gown

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMMS



Disposable Surgical Gown SMMS, Non-Sterile



Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products / EN 13795-1:2019; AAMI PB70: 2012

AAMI liquid barrier level: Level 2

Feature:

Disposable surgical isolation gown-standard performance

SMMS Fabric

Liquid penetration

Scope of Application: It is used for general isolation in clinics, wards, and laboratories of medical institutions.

It is used for routine protection, laboratory protection, work and production resumption protection, and protection during public place disinfection

Structural Composition: The product is mainly made of non-woven cloth through cutting and sewing.

It is provided in a non-sterile condition and only for single use.

Clothing Features and Overview:

The isolation gown is one-piece, coming in both medium and long styles, as shown in the figure.

Rationally designed for free movement, it is easy to put on and take off. The sleeve cuffs have an elastic closure.

It is dry, clean, mildew-free, and free of adhesions, cracks and cavities.

Humanized design: The waist ties with a knot, which is convenient for operating personnel of similar height but different waist size.

Characteristics

Comfortable fabric: New materials, soft, light, breathable.

Waist tie: Use of a knot-tie design allows for flexible adjustment, enhancing the clothing's personalized fit.

Cuff lock-stitching: Elastic cuff design ensures both heightened protection and free movement.

Application:

Our surgical gown is to be worn by medical staff during surgery.

The fabric feature is blue SMMS fabric.

The exam gowns are durable, soft, light, comfortable, abrasion and static electricity resistant.

The exam gowns are air and vapor permeable, while resisting water fluid and aerosol.

The gown is made from 100% high density non woven, containing no additives, adhesives.

The apron features a loop neck closure for a custom fit without gaps.

Cuffs at the wrists for added protection.

All-in-one, front opening

Long sleeved gown with open back, under knee length

Used in surgical/operating room

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).